

ResolutionMD™ Web 2.5

Feature List

SUPPORTED IMAGING MODALITIES

- CT
- MR
- SC
- Enhanced MR
- XA
- CR/DR
- PET
- NM

2D

- Pan
- Zoom
- Window Width
- Window Level
- WL Presets (With Keyboard Shortcuts)
- Image Flip (Horizontal/Vertical)
- Image Rotate (Clockwise/Counterclockwise)
- Text Annotation
- Freehand Annotation
- Point Measurement
- Linear Measurement
- ROI Measurement
- Angle Measurement
- Measurement Calipers
- Magnifying Glass
- Orientation Labels
- Onscreen Image Details (Metadata)
- Reset To Original View Settings
- Network Performance Display
- Screenshot (Send To PACS)
- Comparison Viewer
- Link Views (Window/Level, Pan, Zoom, Etc.)
- Keyboard Shortcut Guide

MIP/MPR

- **All features of 2D except image flip/rotate**
- Oblique MPR
- Orthogonal MPR Reference Lines/Controls
- Slab MPR (Up To 50 mm)
- MIP, minIP, Average Slab Rendering
- Cine
- Reformatted Series Export (With Adjustable Extents And Interval)
- Curved MPR (Including Slab Rendering Modes)
- Volume Measurement
- Orientation Widget

3D

- Pan
- Zoom
- Window Width
- Window Level
- Customizable WL Presets
- Text Annotation
- Freehand Annotation
- Linear Measurement
- Onscreen Image Details (Metadata)
- Reset To Original View Settings
- Network Performance Display
- Screenshot (Send To PACS)
- Keyboard Shortcut Guide
- Orientation Widget
- GPU-based Volume Rendering
- Clipping Planes
- Lens Tool
- Scalpel Tool
- Customizable VR Presets
- Orientation Preset Hotkeys
- Interactive Orientation Widget
- Bone Removal
- Material Editor
- Rotational Batch Capture